

## Serine racemase

Catalog No: #AP88304



Package Size: #AP88304-1 50ug #AP88304-2 100ug #AP88304-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Serine racemase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	D-serine ammonia-lyase,D-serine dehydratase,L-serine ammonia-lyase,L-serine dehydratase
Accession No.	Uniprot ID:Q9QZX7
Uniprot	Q9QZX7
GeneID	27364;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Srr

Note: This product is for in vitro research use only