

# Transmembrane protease serine 4

Catalog No: #AP88309

Package Size: #AP88309-1 50ug #AP88309-2 100ug #AP88309-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |   |
|--------------------|---|
| Product Name       | Transmembrane protease serine 4   |
| Brief Description  | Recombinant Protein   |
| Host Species       | E.coli  |
| Purification       | Greater than 90% by SDS-PAGE  |
| Species Reactivity | Hu  |
| Other Names        | Channel-activating protease 2,CAPH2,Membrane-type serine protease 2,MT-SP2,TMPRSS3,UNQ776/PRO1570 |
| Accession No.      | Uniprot ID:Q9NRS4   |
| Uniprot            | Q9NRS4  |
| GeneID             | 56649;  |
| Tag Info           | His   |
| Storage            | Store at -20°C. (Avoid repeated freezing and thawing.)  |

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:TMPRSS4

Note: This product is for in vitro research use only