

NOLA1 Antibody

Catalog No: #44503

Package Size: #44503-1 50ul #44503-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	NOLA1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	NOLA1 Antibody detects endogenous levels of total NOLA1
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human NOLA1
Target Name	NOLA1
Other Names	GAR1, NOLA1, SnoRNP protein GAR1
Accession No.	Swiss-Prot#: Q9NY12NCBI Gene ID:54433
Uniprot	Q9NY12
GeneID	54433;
Calculated MW	22 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Note: This product is for in vitro research use only