

Regenerating islet-derived protein 4

Catalog No: #AP88314



Package Size: #AP88314-1 50ug #AP88314-2 100ug #AP88314-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Regenerating islet-derived protein 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	REG-4
Accession No.	Uniprot ID:Q9D8G5
Uniprot	Q9D8G5
GeneID	67709;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Reg4

Note: This product is for in vitro research use only