

Relaxin-3

Catalog No: #AP88315



Package Size: #AP88315-1 50ug #AP88315-2 100ug #AP88315-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Relaxin-3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Insulin-like peptide INSL7, Insulin-like peptide 7, Prorelaxin H3, INSL7, RXN3, ZINS4, UNQ6188/PRO20213
Accession No.	Uniprot ID:Q8WXF3
Uniprot	Q8WXF3
GeneID	117579;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: RLN3

Note: This product is for in vitro research use only