

Tumor necrosis factor receptor superfamily member 5

Catalog No: #AP88318



Package Size: #AP88318-1 50ug #AP88318-2 100ug #AP88318-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Tumor necrosis factor receptor superfamily member 5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	B-cell surface antigen CD40,CD40L receptor,CD40,TNFRSF5
Accession No.	Uniprot ID:Q28203
Uniprot	Q28203
GeneID	286849;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CD40

Note: This product is for in vitro research use only