

NAD-dependent ADP-ribosyltransferase sirtuin-4

Catalog No: #AP88321



Package Size: #AP88321-1 50ug #AP88321-2 100ug #AP88321-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	NAD-dependent ADP-ribosyltransferase sirtuin-4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	NAD-dependent ADP-ribosyltransferase sirtuin-4,NAD-dependent protein deacetylase sirtuin-4,Regulatory protein SIR2 homolog 4,SIR2-like protein 4,SIR2L4
Accession No.	Uniprot ID:Q9Y6E7
Uniprot	Q9Y6E7
GeneID	23409;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SIRT4

Note: This product is for in vitro research use only