

Metalloproteinase inhibitor 3

Catalog No: #AP88323



Package Size: #AP88323-1 50ug #AP88323-2 100ug #AP88323-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Metalloproteinase inhibitor 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	21 kDa protein of extracellular matrix, Tissue inhibitor of metalloproteinases 3, TIMP-3, IMP-3
Accession No.	Uniprot ID:P26652
Uniprot	P26652
GeneID	396483;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: TIMP3

Note: This product is for in vitro research use only