

Telomerase protein component 1

Catalog No: #AP88328



Package Size: #AP88328-1 50ug #AP88328-2 100ug #AP88328-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Telomerase protein component 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Telomerase-associated protein 1, Telomerase protein 1, p240, p80 telomerase homolog, Tp1
Accession No.	Uniprot ID:P97499
Uniprot	P97499
GeneID	21745;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Tep1

Note: This product is for in vitro research use only