

Tissue factor

Catalog No: #AP88329



Package Size: #AP88329-1 50ug #AP88329-2 100ug #AP88329-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Tissue factor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	TF,Coagulation factor III,CD142
Accession No.	Uniprot ID:P30931
Uniprot	P30931
GeneID	280686;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:F3

Note: This product is for in vitro research use only