

Mothers against decapentaplegic homolog 6

Catalog No: #AP88332



Package Size: #AP88332-1 50ug #AP88332-2 100ug #AP88332-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mothers against decapentaplegic homolog 6
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	MAD homolog 6, Mad homolog 7, SMAD family member 6, SMAD 6, Madh6, Madh7, Msmad6
Accession No.	Uniprot ID:O35182
Uniprot	O35182
GeneID	17130;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Smad6

Note: This product is for in vitro research use only