

# Tryptophan 5-hydroxylase 1

Catalog No: #AP88335



Package Size: #AP88335-1 50ug #AP88335-2 100ug #AP88335-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Tryptophan 5-hydroxylase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Tryptophan 5-monoxygenase 1, Tph
Accession No.	Uniprot ID:P09810
Uniprot	P09810
GeneID	24848;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Tph1

Note: This product is for in vitro research use only