

Transient receptor potential cation channel subfamily V member 1

Catalog No: #AP88336

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP88336-1 50ug #AP88336-2 100ug #AP88336-3 1mg

Description

Product Name	Transient receptor potential cation channel subfamily V member 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	TrpV1,Osm-9-like TRP channel 1,OTRPC1,Vanilloid receptor 1
Accession No.	Uniprot ID:Q704Y3
Uniprot	Q704Y3
GeneID	193034;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Trpv1

Note: This product is for in vitro research use only