

Sorting nexin-9

Catalog No: #AP88339



Package Size: #AP88339-1 50ug #AP88339-2 100ug #AP88339-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Sorting nexin-9
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Accession No.	Uniprot ID:Q91VH2
Uniprot	Q91VH2
GeneID	66616;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Snx9

Note: This product is for in vitro research use only