

Tumor necrosis factor receptor superfamily member 10B

Catalog No: #AP88340

Package Size: #AP88340-1 50ug #AP88340-2 100ug #AP88340-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Tumor necrosis factor receptor superfamily member 10B
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Death receptor 5, TNF-related apoptosis-inducing ligand receptor 2, TRAIL receptor 2, CD262, DR5, KILLER, TRAILR2, TRICK2, ZTNFR9, UNQ160/PRO186
Accession No.	Uniprot ID:O14763
Uniprot	O14763
GeneID	8795;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: TNFRSF10B

Note: This product is for in vitro research use only