

## RISC-loading complex subunit TARBP2

Catalog No: #AP88341



Package Size: #AP88341-1 50ug #AP88341-2 100ug #AP88341-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	RISC-loading complex subunit TARBP2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Protamine-1 RNA-binding protein,PRM-1 RNA-binding protein,TAR RNA-binding protein 2,Prbp
Accession No.	Uniprot ID:P97473
Uniprot	P97473
GeneID	21357;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Tarbp2

Note: This product is for in vitro research use only