

Serine/threonine-protein kinase TAO1

Catalog No: #AP88342



Package Size: #AP88342-1 50ug #AP88342-2 100ug #AP88342-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Serine/threonine-protein kinase TAO1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Thousand and one amino acid protein 1,Kiaa1361
Accession No.	Uniprot ID:Q5F2E8
Uniprot	Q5F2E8
GeneID	216965;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Taok1

Note: This product is for in vitro research use only