

T-cell leukemia/lymphoma protein 1A

Catalog No: #AP88343



Package Size: #AP88343-1 50ug #AP88343-2 100ug #AP88343-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	T-cell leukemia/lymphoma protein 1A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Oncogene TCL-1,Oncogene TCL1,Protein p14 TCL1,Tcl1
Accession No.	Uniprot ID:P56280
Uniprot	P56280
GeneID	21432;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Tcl1a

Note: This product is for in vitro research use only