

TAR DNA-binding protein 43

Catalog No: #AP88344



Package Size: #AP88344-1 50ug #AP88344-2 100ug #AP88344-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	TAR DNA-binding protein 43
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	TDP-43,Tdp43
Accession No.	Uniprot ID:Q921F2
Uniprot	Q921F2
GeneID	230908;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Tardbp

Note: This product is for in vitro research use only