

Switch-associated protein 70

Catalog No: #AP88349



Package Size: #AP88349-1 50ug #AP88349-2 100ug #AP88349-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Switch-associated protein 70
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	SWAP-70,RCJMB04_1e1,RCJMB04_8b3
Accession No.	Uniprot ID:Q5F4B2
Uniprot	Q5F4B2
GeneID	423044;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SWAP70

Note: This product is for in vitro research use only