

Tumor necrosis factor alpha-induced protein 3

Catalog No: #AP88350



Package Size: #AP88350-1 50ug #AP88350-2 100ug #AP88350-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Tumor necrosis factor alpha-induced protein 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	TNF alpha-induced protein 3,Putative DNA-binding protein A20,Zinc finger protein A20,Tnfip3
Accession No.	Uniprot ID:Q60769
Uniprot	Q60769
GeneID	21929;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Tnfaip3

Note: This product is for in vitro research use only