

Tumor necrosis factor ligand superfamily member 13

Catalog No: #AP88351



Package Size: #AP88351-1 50ug #AP88351-2 100ug #AP88351-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Tumor necrosis factor ligand superfamily member 13
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	A proliferation-inducing ligand,APRIL,TNF- and APOL-related leukocyte expressed ligand 2,TALL-2,TNF-related death ligand 1,TRDL-1,CD256,APRIL,TALL2,ZTNF2,UNQ383/PRO715
Accession No.	Uniprot ID:O75888
Uniprot	O75888
GeneID	8741;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:TNFSF13

Note: This product is for in vitro research use only