

## Tumor necrosis factor ligand superfamily member 6

Catalog No: #AP88352



Package Size: #AP88352-1 50ug #AP88352-2 100ug #AP88352-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Tumor necrosis factor ligand superfamily member 6
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	CD95 ligand,CD95-L,Fas antigen ligand,Fas ligand,CD178,Apt1Lg1,Cd95l,Fasl,Tnfsf6
Accession No.	Uniprot ID:P36940
Uniprot	P36940
GeneID	25385;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Faslg

Note: This product is for in vitro research use only