

Tenomodulin

Catalog No: #AP88353



Package Size: #AP88353-1 50ug #AP88353-2 100ug #AP88353-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Tenomodulin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	TeM, Chondromodulin-1-like protein, ChM1L, Chondromodulin-I-like protein, Myodulin, Tendin, CHM1L, UNQ771/PRO1565
Accession No.	Uniprot ID:Q9H2S6
Uniprot	Q9H2S6
GeneID	64102;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: TNMD

Note: This product is for in vitro research use only