

Transferrin receptor protein 1

Catalog No: #AP88354



Package Size: #AP88354-1 50ug #AP88354-2 100ug #AP88354-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Transferrin receptor protein 1 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Ms |
| Other Names | TR,CD71,Trfr |
| Accession No. | Uniprot ID:Q62351 |
| Uniprot | Q62351 |
| GeneID | 22042; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Tfrc

Note: This product is for in vitro research use only