

# Transforming growth factor beta-3

Catalog No: #AP88355



Package Size: #AP88355-1 50ug #AP88355-2 100ug #AP88355-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Transforming growth factor beta-3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	TGF-beta-3,Tgf-b3
Accession No.	Uniprot ID:Q07258
Uniprot	Q07258
GeneID	25717;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Tgfb3

Note: This product is for in vitro research use only