

## Tumor necrosis factor ligand superfamily member 11

Catalog No: #AP88360



Package Size: #AP88360-1 50ug #AP88360-2 100ug #AP88360-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Tumor necrosis factor ligand superfamily member 11
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Osteoclast differentiation factor,ODF,Osteoprotegerin ligand,OPGL,Receptor activator of nuclear factor kappa-B ligand,RANKL,TNF-related activation-induced cytokine,TRANCE,CD254,Opgl,Rankl,Trance
Accession No.	Uniprot ID:Q9ESE2
Uniprot	Q9ESE2
GeneID	117516;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Tnfsf11

Note: This product is for in vitro research use only