

Tumor necrosis factor receptor superfamily member 1B

Catalog No: #AP88361

Package Size: #AP88361-1 50ug #AP88361-2 100ug #AP88361-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Tumor necrosis factor receptor superfamily member 1B
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Tumor necrosis factor receptor 2, TNF-R2, Tumor necrosis factor receptor type II, TNF-RII, p75, p80 TNF-alpha receptor, CD120b, Tnfr-2, Tnfr2
Accession No.	Uniprot ID:P25119
Uniprot	P25119
GeneID	21938;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Tnfrsf1b

Note: This product is for in vitro research use only