

Transcription factor 21

Catalog No: #AP88365



Package Size: #AP88365-1 50ug #AP88365-2 100ug #AP88365-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Transcription factor 21
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	TCF-21,Capsulin,Class A basic helix-loop-helix protein 23,bHLHa23,Epicardin,Podocyte-expressed 1,Pod-1,BHLHA23,POD1
Accession No.	Uniprot ID:O43680
Uniprot	O43680
GeneID	6943;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:TCF21

Note: This product is for in vitro research use only