

Transcobalamin-2

Catalog No: #AP88366



Package Size: #AP88366-1 50ug #AP88366-2 100ug #AP88366-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Transcobalamin-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	TC-2,Transcobalamin II,TC II
Accession No.	Uniprot ID:O88968
Uniprot	O88968
GeneID	21452;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Tcn2

Note: This product is for in vitro research use only