

Transforming growth factor beta-1

Catalog No: #AP88367



Package Size: #AP88367-1 50ug #AP88367-2 100ug #AP88367-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Transforming growth factor beta-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	TGF-beta-1
Accession No.	Uniprot ID:P09531
Uniprot	P09531
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:TGFB1

Note: This product is for in vitro research use only