

Cytochrome P450 11B2 Antibody

Catalog No: #44509

Package Size: #44509-1 50ul #44509-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Cytochrome P450 11B2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	Cytochrome P450 11B2 Antibody detects endogenous levels of total Cytochrome P450 11B2
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human Cytochrome P450 11B2
Target Name	Cytochrome P450 11B2
Other Names	CYP11B2, ALDOS, Aldosterone synthase, CYP11BL, CYPXIB2, Cytochrome P-450Aldo, Cytochrome P-450C18, P-450C18, p450aldo, Steroid 18-hydroxylase, Steroid 11-beta-monooxygenase, CYP11B, p450C18, Steroid 11-beta/18-hydroxylase
Accession No.	Swiss-Prot#: P19099NCBI Gene ID:1585
Uniprot	P19099
GeneID	1585;
Calculated MW	58 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western Blot analysis of various cells using CYP11B1/2 Polyclonal Antibody

Note: This product is for in vitro research use only