

Zinc finger protein SNAI2

Catalog No: #AP88372



Package Size: #AP88372-1 50ug #AP88372-2 100ug #AP88372-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Zinc finger protein SNAI2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Neural crest transcription factor Slug, Protein snail homolog 2, SLUG
Accession No.	Uniprot ID:Q3MHQ4
Uniprot	Q3MHQ4
GeneID	520631;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: SNAI2

Note: This product is for in vitro research use only