

Sphingosine kinase 1

Catalog No: #AP88377



Package Size: #AP88377-1 50ug #AP88377-2 100ug #AP88377-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Sphingosine kinase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	SK 1,SPHK,SPK
Accession No.	Uniprot ID:Q9NYA1
Uniprot	Q9NYA1
GeneID	8877;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SPHK1

Note: This product is for in vitro research use only