

Serine/threonine-protein kinase TAO2

Catalog No: #AP88379



Package Size: #AP88379-1 50ug #AP88379-2 100ug #AP88379-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Serine/threonine-protein kinase TAO2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Kinase from chicken homolog C,hKFC-C,Prostate-derived sterile 20-like kinase 1,PSK-1,Thousand and one amino acid protein kinase 2,KIAA0881,MAP3K17,PSK,PSK1,UNQ2971/PRO7431
Accession No.	Uniprot ID:Q9UL54
Uniprot	Q9UL54
GeneID	9344;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:TAOK2

Note: This product is for in vitro research use only