

TGF-beta-activated kinase 1 and MAP3K7-binding protein 2

Catalog No: #AP88380

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP88380-1 50ug #AP88380-2 100ug #AP88380-3 1mg

Description

Product Name	TGF-beta-activated kinase 1 and MAP3K7-binding protein 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Mitogen-activated protein kinase kinase kinase 7-interacting protein 2, TAK1-binding protein 2, TAB-2, TGF-beta-activated kinase 1-binding protein 2, KIAA0733, MAP3K7IP2
Accession No.	Uniprot ID: Q9NYJ8
Uniprot	Q9NYJ8
GeneID	23118;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: TAB2

Note: This product is for in vitro research use only