

Tubulin alpha-1A chain

Catalog No: #AP88381



Package Size: #AP88381-1 50ug #AP88381-2 100ug #AP88381-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Tubulin alpha-1A chain
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Alpha-tubulin 3,Tubulin B-alpha-1,Tubulin alpha-3 chain,TUBA3
Accession No.	Uniprot ID:Q71U36
Uniprot	Q71U36
GeneID	7846;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:TUBA1A

Note: This product is for in vitro research use only