

Sorbin and SH3 domain-containing protein 1

Catalog No: #AP88391



Package Size: #AP88391-1 50ug #AP88391-2 100ug #AP88391-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Sorbin and SH3 domain-containing protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Ponsin,SH3 domain protein 5,SH3P12,c-Cbl-associated protein,CAP,Kiaa1296,Sh3d5
Accession No.	Uniprot ID:Q62417
Uniprot	Q62417
GeneID	20411;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Sorbs1

Note: This product is for in vitro research use only