

Retinol-binding protein 4

Catalog No: #AP88394



Package Size: #AP88394-1 50ug #AP88394-2 100ug #AP88394-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Retinol-binding protein 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Plasma retinol-binding protein,PRBP
Accession No.	Uniprot ID:P04916
Uniprot	P04916
GeneID	25703;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Rbp4

Note: This product is for in vitro research use only