

# Lymphotoxin-beta

Catalog No: #AP88395



Package Size: #AP88395-1 50ug #AP88395-2 100ug #AP88395-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Lymphotoxin-beta
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	LT-beta,Tumor necrosis factor C,TNF-C,Tumor necrosis factor ligand superfamily member 3,Tnfc,Tnfsf3
Accession No.	Uniprot ID:P41155
Uniprot	P41155
GeneID	16994;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Ltb

Note: This product is for in vitro research use only