

Sperm-specific antigen 2 homolog

Catalog No: #AP88396



Package Size: #AP88396-1 50ug #AP88396-2 100ug #AP88396-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Sperm-specific antigen 2 homolog
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Ki-ras-induced actin-interacting protein,Kiaa1927,Krap
Accession No.	Uniprot ID:Q922B9
Uniprot	Q922B9
GeneID	70599;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ssfa2

Note: This product is for in vitro research use only