

## General transcription factor IIH subunit 2

Catalog No: #AP88398



Package Size: #AP88398-1 50ug #AP88398-2 100ug #AP88398-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	General transcription factor IIH subunit 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Basic transcription factor 2 44 kDa subunit,BTF2 p44,General transcription factor IIH polypeptide 2,TFIIH basal transcription factor complex p44 subunit,BTF2P44
Accession No.	Uniprot ID:Q13888
Uniprot	Q13888
GeneID	2966;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:GTF2H2

Note: This product is for in vitro research use only