

Mothers against decapentaplegic homolog 5

Catalog No: #AP88403



Package Size: #AP88403-1 50ug #AP88403-2 100ug #AP88403-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mothers against decapentaplegic homolog 5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	MAD homolog 5, JV5-1, SMAD family member 5, SMAD 5, MADH5
Accession No.	Uniprot ID:Q99717
Uniprot	Q99717
GeneID	4090;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: SMAD5

Note: This product is for in vitro research use only