

Thyroid hormone receptor beta

Catalog No: #AP88410



Package Size: #AP88410-1 50ug #AP88410-2 100ug #AP88410-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Thyroid hormone receptor beta
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Nuclear receptor subfamily 1 group A member 2,c-erbA-2,c-erbA-beta,ERBA2,NR1A2,THR1
Accession No.	Uniprot ID:P10828
Uniprot	P10828
GeneID	7068;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:THRB

Note: This product is for in vitro research use only