

Serine--tRNA ligase, cytoplasmic

Catalog No: #AP88412



Package Size: #AP88412-1 50ug #AP88412-2 100ug #AP88412-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Serine--tRNA ligase, cytoplasmic
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Seryl-tRNA synthetase, SerRS, Seryl-tRNA(Ser/Sec) synthetase, Sars1, Sers
Accession No.	Uniprot ID:P26638
Uniprot	P26638
GeneID	20226;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Sars

Note: This product is for in vitro research use only