

Transmembrane 4 L6 family member 20

Catalog No: #AP88413



Package Size: #AP88413-1 50ug #AP88413-2 100ug #AP88413-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Transmembrane 4 L6 family member 20
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	UNQ518/PRO994
Accession No.	Uniprot ID:Q53R12
Uniprot	Q53R12
GeneID	79853;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: TM4SF20

Note: This product is for in vitro research use only