

Spectrin alpha chain, brain

Catalog No: #AP88414



Package Size: #AP88414-1 50ug #AP88414-2 100ug #AP88414-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Spectrin alpha chain, brain
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Alpha-II spectrin,Fodrin alpha chain,Spectrin, non-erythroid alpha chain,Spna2,Spta2
Accession No.	Uniprot ID:P16546
Uniprot	P16546
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Sptan1

Note: This product is for in vitro research use only