

Tumor necrosis factor receptor superfamily member 1A

Catalog No: #AP88419

Package Size: #AP88419-1 50ug #AP88419-2 100ug #AP88419-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Tumor necrosis factor receptor superfamily member 1A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Tumor necrosis factor receptor 1, TNF-R1, Tumor necrosis factor receptor type I, TNF-R1, p55, p60, CD120a, Tnfr-1, Tnfr1
Accession No.	Uniprot ID:P22934
Uniprot	P22934
GeneID	25625;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Tnfrsf1a

Note: This product is for in vitro research use only