

## Sclerostin domain-containing protein 1

Catalog No: #AP88422



Package Size: #AP88422-1 50ug #AP88422-2 100ug #AP88422-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Sclerostin domain-containing protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Ectodermal BMP inhibitor,Ectodin,Uterine sensitization-associated gene 1 protein,USAG-1,USAG1,CDA019
Accession No.	Uniprot ID:Q6X4U4
Uniprot	Q6X4U4
GeneID	25928;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:SOSTDC1

Note: This product is for in vitro research use only