

DNA nucleotidylexotransferase

Catalog No: #AP88424



Package Size: #AP88424-1 50ug #AP88424-2 100ug #AP88424-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	DNA nucleotidylexotransferase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	Terminal addition enzyme, Terminal deoxynucleotidyltransferase, Terminal transferase, TDT
Accession No.	Uniprot ID:P36195
Uniprot	P36195
GeneID	396351;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: DNNT

Note: This product is for in vitro research use only